



Figure 35. Plasmid mixtures purified from the *Hind* III genomic DNA sub-library. Plasmids were digested by *Hind* III and analyzed on a 0.8% agarose gel. Lanes 1 and 32: Lambda *Sy* I DNA standard; lanes 2 to 16: plasmid mixtures purified from pooled row tubes #1 to # 15; lanes 17 to 31: Plasmid mixtures purified from pooled column tubes # A to # O. Lane 7: pooled row tube# 6 was selected for the further screening.

[previous figure][next figure]