

## *Appendix 2*

### **Storage of Pythiaceae on Hemp Seed**

---

Use 19 x 65 mm glass vials with screw-on cap with rubber-lined closure or similar glass vial and fill with approximately 5 ml dH<sub>2</sub>O. Place two hemp seeds in each vial; cap loosely and autoclave. When cool tighten cap. Store at 4°C. Place 3 to 5 mm agar plugs of a pure culture in vial. Place cap on tightly and after active mycelial growth appears (1 to 2 days), seal with parafilm. Store at room temperature. Recover culture by placing colonized agar plugs or hemp seeds on an appropriate culture medium.