

PE	1	2
	<b>Single Positive for PE</b>	<b>Double Positive</b>
	3	4
	<b>Double Negative</b>	<b>Single Positive for FITC</b>
	<b>FITC</b>	

**Appendix Figure 2. Illustration of how double-stained cell samples are analyzed.** Cells are analyzed flow cytometrically for both FITC and PE fluorescence. Cells staining positively for both stains appear in the second quadrant. Single positive cells appear in either the first or fourth quadrants, and double-negative cells appear in the third quadrant.